

GOAT ANTI-MUNC18C (AA514-527) ANTIBODY

SKU: EB12335



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	STXBP3 syntaxin binding protein 3 MUNC18-3 MUNC18C PSP UNC-18C platelet Sec1 protein protein unc-18 homolog 3 protein unc-18 homolog C syntaxin 4 binding protein syntaxin-binding protein 3 unc18-3
Accession ID	NP_009200.2
Blocking Peptide	EBP12335
Immunogen	Peptide with sequence C-RQKPRANYLEDKRN, from the internal region of the protein sequence according to NP_009200.2.
Peptide Sequence	C-RQKPRANYLEDKRN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Pig
Reactive Species	Human
Human Gene ID	6814
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 60kDa band observed in Human Platelets lysates (calculated MW of 67.8kDa according to NP_009200.2). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

