

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-MUPCDH ANTIBODY

SKU: EB08372

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

MUPCDH|mucin-like protocadherin|FLJ20219|MU-PCDH|MUCDHL|mu-protocadherin|mucin and cadherin-like|mucin **Alias**

and cadherin-like protein **Names**

Accession

NP 068743.2; NP 112554.2

Blocking

EBP08372 **Peptide**

Peptide with sequence C-PYDAPGGDDSY, from the C Terminus of the protein sequence according to NP 068743.2; Immunogen

NP 112554.2.

Product Comments

This antibody is expected to recognise both reported isoforms (NP_068743.2; NP_112554.2).

Peptide

C-PYDAPGGDDSY

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human Species

Human

53841

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:8000.

Limit

Western Blot

Preliminary experiments in Human Colon, Ileum and Kidney lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any

results been reported with other antibodies/lysates?

Application

Pep-ELISA **Type**





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600