

GOAT ANTI-MUPCDH ANTIBODY

SKU: EB08372



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	MUPCDH mucin-like protocadherin FLJ20219 MU-PCDH MUCDHL mu-protocadherin mucin and cadherin-like mucin and cadherin-like protein
Accession ID	NP_068743.2; NP_112554.2
Blocking Peptide	EBP08372
Immunogen	Peptide with sequence C-PYDAPGGDDSY, from the C Terminus of the protein sequence according to NP_068743.2; NP_112554.2.
Product Comments	This antibody is expected to recognise both reported isoforms (NP_068743.2; NP_112554.2).
Peptide Sequence	C-PYDAPGGDDSY
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	53841
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Preliminary experiments in Human Colon, Ileum and Kidney lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)