


GOAT ANTI-MURF1 / TRIM63 (N TERM) ANTIBODY

SKU: EB05879



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage	
Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	muscle specific ring finger protein 1 tripartite motif-containing 63 striated muscle RING zinc finger protein muscle specific ring finger protein 2 iris ring finger protein FLJ32380 MURF2 SMRZ IRF ring finger protein
Names	28 RNF28 MURF-1 MURF1 TRIM63
Usage Summary	Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A549 cells. Flow Cytometry: Flow cytometric analysis of A549 cells. Recommended concentration: 10ug/ml.
Accession ID	NP_115977.2
Blocking Peptide	EBP05879
Immunogen	Peptide with sequence DYKSSLIQDGNPM-C, from the N Terminus of the protein sequence according to NP_115977.2.

Peptide Sequence	DYKSSLIQDGNPM-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Cow
Reactive Species	Human
Human Gene ID	84676
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png
IHC Results	In paraffin embedded Human Skeletal Muscle shows strong staining in rare nuclei, and weak staining of muscle fibres in longitudinal section. Recommended concentration: 1.25µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, IHC, IF, FC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

