




Telephone: (650) 697-3600



GOAT ANTI-MURF1 / TRIM63 (N TERM) ANTIBODY

SKU: EB05879



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / muscle specific ring finger protein 1|tripartite motif-containing 63|striated muscle RING zinc finger protein|muscle

Alias specific ring finger protein 2|iris ring finger protein|FLJ32380|MURF2|SMRZ|IRF|ring finger protein

Names 28|RNF28|MURF-1|MURF1|TRIM63

Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A549 cells. Usage Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of A549

Summary

cells. Recommended concentration: 10ug/ml.

Accession NP 115977.2

Blocking EBP05879 **Peptide**

Immunogen Peptide with sequence DYKSSLIQDGNPM-C, from the N Terminus of the protein sequence according to NP 115977.2.









Telephone: (650) 697-3600

Peptide DYKSSLIQDGNPM-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Refrigerated Instructions

Predicted

Human, Cow **Species**

Reactive

Human **Species**

Human

84676

Gene ID Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png Grade

In paraffin embedded Human Skeletal Muscle shows strong staining in rare nuclei, and weak staining of muscle **IHC Results**

fibres in longitudinal section. Recommended concentration: 1.25µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:32000.

Limit

Application

Type

Pep-ELISA, IHC, IF, FC

GALLERY IMAGES

