

GOAT ANTI-MYLIP ANTIBODY

SKU: EB09591

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	cellular modulator of immune recognition (c-MIR) band 4.1 superfamily member
Names	BZF1 OTTHUMP00000018026 MIR myosin regulatory light chain interacting protein MYLIP
Accession ID	NP_037394.2
Blocking Peptide	EBP09591
Immunogen	Peptide with sequence QNTAKYNYEELCAKE, from the internal region of the protein sequence according to NP_037394.2.
Peptide Sequence	QNTAKYNYEELCAKE
Shipping Instructions	Refrigerated
Predicted Species	Human; Dog, Cow, Pig
Reactive Species	Human
Human Gene ID	29116
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 49kDa band observed in lysates of cell line HeLa (calculated MW of 49.9kDa according to NP_037394.2). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

