

GOAT ANTI-MYO18B ANTIBODY

SKU: EB05849



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	FLJ41942 DKFZp779C1668 DKFZp434C2422 CTA-109D1.1 myosin 18B bK125H2.1 BK125H2.1 myosin XVIIIIB MYO18B
Accession ID	NP_115997.5
Blocking Peptide	EBP05849
Immunogen	Peptide with sequence C-DDVASIMKKYLQK, from the C Terminus of the protein sequence according to NP_115997.5.
Peptide Sequence	C-DDVASIMKKYLQK
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	84700
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Preliminary experiments gave no signal but low background in Human Heart and A431 lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA