

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-MYOSIN IXB ANTIBODY

SKU: EB07266



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa 20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

celiac disease 4|myosin|MYR5|CELIAC4|myosin IXB|MYO9B

Names

Accession

NP_004136.2; NP_001123537.1

Blocking

EBP07266 **Peptide**

Peptide with sequence C-KNDFEPVKQSKAQ, from the internal region of the protein sequence according to **Immunogen**

NP_004136.2; NP_001123537.1.

Product Comments

This antibody is expected to recognise both reported isoforms (NP 004136.2; NP 001123537.1).

Peptide

Sequence

C-KNDFEPVKQSKAQ

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human

Species

Human 4650

Gene ID **Product**

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:2000.

Limit

Western Approx 200kDa band observed in Human Liver lysates (calculated MW of 243kDa according to NP_004136.2).

Recommended concentration: 2-6µg/ml. Blot

Application

Pep-ELISA, WB Type

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa 20kDa

15kDa