

GOAT ANTI-MYTH15 ANTIBODY

SKU: EB09748



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	myosin, heavy polypeptide 15 myosin, heavy chain 15 MYH15
Accession ID	NP_055796.1
Blocking Peptide	EBP09748
Immunogen	Peptide with sequence TADGESLSIKEDKIQ-C, from the N Terminus of the protein sequence according to NP_055796.1.
Peptide Sequence	TADGESLSIKEDKIQ-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	22989
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)