Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-MYTH15 ANTIBODY

SKU: EB09748

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias myosin, heavy polypeptide 15|myosin, heavy chain 15|MYH15

Names

Accession NP 055796.1 ID

Blocking

EBP09748 **Peptide**

Immunogen NP_055796.1. Peptide with sequence TADGESLSIKEDKIQ-C, from the N Terminus of the protein sequence according to

Peptide

TADGESLSIKEDKIQ-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

22989 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:4000.

Limit

Application Pep-ELISA **Type**