

GOAT ANTI-NAPSA / NAPSIN A ANTIBODY

SKU: EB05032



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	napsin A napsin A aspartic peptidase SNAPA NAPA NAP1 Kdap KAP pronapsin A
Accession ID	NP_004842.1
Blocking Peptide	EBP05032
Immunogen	Peptide with sequence C-LGWGETAQAQFPG, from the C Terminus of the protein sequence according to NP_004842.1.
Peptide Sequence	C-LGWGETAQAQFPG
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	9476
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western blot: Approx 40 kDa band seen in 293 Lysate. Recommended for use at 1-3 µg/ml.
Application Type	Pep-ELISA, WB

