

GOAT ANTI-NAT2 (AA177-188) ANTIBODY

SKU: EB07244



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	arylamide acetylase 2 (N-acetyltransferase 2, isoniazid inactivation) arylamide acetylase 2 Arylamine N-acetyltransferase-2 AAC2 HGNC:7646 N-acetyltransferase 2 (arylamine N-acetyltransferase) NAT2
Accession ID	NP_000006.1
Blocking Peptide	EBP07244
Immunogen	Peptide with sequence C-NSHLLPKKKHQK, from the internal region of the protein sequence according to NP_000006.1.
Peptide Sequence	C-NSHLLPKKKHQK
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	10
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Preliminary experiments in human colon, duodenum and liver lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA