



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-NBR2 ANTIBODY

SKU: EB05162

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias neighbor of BRCA1 gene 2|NBR2

Names

Accession NP 005812

Immunogen Peptide with sequence C-TQGFLNPLEFSA, from the C Terminus of the protein sequence according to NP_005812.

Peptide C-TQGFLNPLEFSA

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Instructions Refrigerated

Predicted

Human

Species

Human

Gene ID

10230

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Western

Western blot: Weak band observed in 293 at approx 30-40kDa with 1 hour incubation of primary. Blot