

# GOAT ANTI-NBR2 ANTIBODY

**SKU:** EB05162



## SPECIFICATIONS

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

**Storage Instructions** Aliquot and store at -20°C. Minimize freezing and thawing.

**Synonym /**

**Alias** neighbor of BRCA1 gene 2|NBR2

**Names**

**Accession ID** NP\_005812

**Immunogen** Peptide with sequence C-TQGFLNPLEFSA, from the C Terminus of the protein sequence according to NP\_005812.

**Peptide Sequence** C-TQGFLNPLEFSA

**Purification Method** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

**Shipping Instructions** Refrigerated

**Predicted Species** Human

**Human Gene ID** 10230

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png)

**Western Blot** Western blot: Weak band observed in 293 at approx 30-40kDa with 1 hour incubation of primary.