

GOAT ANTI-NCBP1 / CBP80 ANTIBODY

SKU: EB07074



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	MGC2087 RP11-546O6.1 nuclear cap binding protein, 80kD nuclear cap binding protein 1, 80kD NCBP CBP80 HGNC:7658 nuclear cap binding protein subunit 1, 80kDa NCBP1
Accession ID	NP_002477.1
Blocking Peptide	EBP07074
Immunogen	Peptide with sequence C-KELYEKDAEMD, from the internal region of the protein sequence according to NP_002477.1.
Peptide Sequence	C-KELYEKDAEMD
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	4686
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Preliminary experiments in lysates of cell lines HEK293, HeLa and HepG2 gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA