

GOAT ANTI-NEDD8 CONJUGATING ENZYME (N TERMINUS) ANTIBODY

SKU: EB10869

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	UBE2F ubiquitin-conjugating enzyme E2F (putative) MGC18120 NCE2 NEDD8 carrier protein UBE2F NEDD8 conjugating enzyme NEDD8 protein ligase UBE2F NEDD8-conjugating enzyme 2 NEDD8-conjugating enzyme UBE2F OTTHUMP00000202755 OTTH
Accession ID	NP_542409.1
Blocking Peptide	EBP10869
Immunogen	Peptide with sequence SKLKRDDGLKGSRT-C, from the N Terminus of the protein sequence according to NP_542409.1.
Peptide Sequence	SKLKRDDGLKGSRT-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Reactive Species	Rat
Human Gene ID	140739
Mouse Gene ID	67921
Rat Gene ID	363284
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 22kDa band observed in Rat Pancreas lysates (calculated MW of 21.1kDa according to Human NP_542409.1 and 18.6kDa according to Rat NP_001008382.1). Recommended concentration: 0.05-0.2µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa