

GOAT ANTI-NEIL1 ANTIBODY

SKU: EB07374



SPECIFICATIONS

Unit Size	100 µg
Synonym /	
Alias	nei endonuclease VIII-like 1 endonuclease VIII hFPG1 NEI1 FLJ22402 nei endonuclease VIII-like 1 (E. coli) NEIL1
Names	
Accession ID	NP_078884.1
Blocking Peptide	EBP07374
Immunogen	Peptide with sequence C-KVKADIPSLEPE, from the C Terminus of the protein sequence according to NP_078884.1.
Peptide Sequence	C-KVKADIPSLEPE
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	79661
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Western Blot: Preliminary experiments gave an approx 75kDa band in Human Liver lysates after 0.01µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 43.7kDa according to NP_078884.1. The 75kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?
Application Type	Pep-ELISA