



Telephone: (650) 697-3600

GOAT ANTI-NESCA ADAPTER PROTEIN ANTIBODY

SKU: EB05085

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

new molecule containing SH3 at the carboxy-terminus|nesca protein|RUN and SH3 domain containing **Alias**

1|NESCA|RUSC1 **Names**

Accession

NP 055143 ID

Peptide with sequence C-EGLVPVGYTSLVL, from the C Terminus of the protein sequence according to **Immunogen**

NP 055143.

Peptide

C-EGLVPVGYTSLVL Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive **Species**

Human

Human

23623

Gene ID Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

Western

Western blot: Approx 50kDa band seen in HeLa and U937 cell lysates. Recommended for use at 1-3 μg/ml. **Blot**