

GOAT ANTI-NESCA ADAPTER PROTEIN ANTIBODY

SKU: EB05085



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	new molecule containing SH3 at the carboxy-terminus nesca protein RUN and SH3 domain containing 1 NESCA RUSC1
Accession ID	NP_055143
Immunogen	Peptide with sequence C-EGLVPVGYTSLVL, from the C Terminus of the protein sequence according to NP_055143.
Peptide Sequence	C-EGLVPVGYTSLVL
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	23623
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 50kDa band seen in HeLa and U937 cell lysates. Recommended for use at 1-3 µg/ml.