

GOAT ANTI-NETRIN-G1 ANTIBODY

SKU: EB12049



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	netrin-G1 netrin G1f laminet-1 laminet 1 axon guidance molecule YLSR571 Lmnt1 netrin G1 NTNG1
Accession ID	NP_001106697.1; NP_001106699.1; NP_055732.2
Blocking Peptide	EBP12049
Immunogen	Peptide with sequence C-SLYGQLDTTKLRD, from the intenal region of the protein sequence according to NP_001106697.1; NP_001106699.1; NP_055732.2.
Product Comments	This antibody is expected to recognize 1, 2, 3.
Peptide Sequence	C-SLYGQLDTTKLRD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	22854
Mouse Gene ID	80883
Rat Gene ID	295382
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)