


GOAT ANTI-NEURTURIN ANTIBODY

SKU: EB07938



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias NTN|neurturin|NRTN

Names

Accession ID NP_004549.1

Blocking Peptide EBP07938

Immunogen Peptide with sequence C-DAHSRYHTVHELARE, from the C Terminus of the protein sequence according to NP_004549.1.

Peptide Sequence C-DAHSRYHTVHELARE

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human
Reactive Species	Human
Human Gene ID	4902
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Colon (Submucosal Plexus). Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 25kDa band observed in Human Heart lysates (calculated MW of 22.4kDa according to NP_004549.1). Recommended concentration: 0.1-0.3µg/ml. Primary incubation was 1 hour.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

