




GOAT ANTI-NEURTURIN ANTIBODY

SKU: EB07938



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias NTN|neurturin|NRTN

Names

ID

Accession

NP_004549.1

Blocking

EBP07938

Peptide

Peptide with sequence C-DAHSRYHTVHELSARE, from the C Terminus of the protein sequence according to Immunogen NP 004549.1.

Peptide

Sequence

C-DAHSRYHTVHELSARE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted

Human

Species Reactive

Human

Species Human

Gene ID

4902

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results Paraffin embedded Human Colon (Submucosal Plexus). Recommended concentration: 5µg/ml.

ELISA

Detection Antibody detection limit dilution 1:64000.

Limit

Approx 25kDa band observed in Human Heart lysates (calculated MW of 22.4kDa according to NP_004549.1).

Blot Recommended concentration: 0.1-0.3μg/ml. Primary incubation was 1 hour.

Application

Western

Type


Pep-ELISA, WB, IHC

GALLERY IMAGES



Telephone: (650) 697-3600





250kDa 150kDa 100kDa 75kDa 50kDa 37kDa 25kDa 20kDa

15kDa