

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-NEURTURIN (MOUSE) ANTIBODY**

**SKU:** EB07939

Telephone: (650) 697-3600





250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa 20kDa

15kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

**Alias** Nrtn|neurturin|NTN

Names

**Accession** 

NP\_032764.1 ID

**Blocking** 

EBP07939 **Peptide** 

Peptide with sequence C-DVHSRYHTLQELSARE, from the C Terminus of the protein sequence according to **Immunogen** 

NP 032764.1.

**Peptide** 

C-DVHSRYHTLQELSARE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Mouse, Rat **Species** 

Reactive

Mouse, Rat **Species** 

Mouse

18188 Gene ID

Rat Gene ID 84423

**Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite\_medium.png$ 

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:64000.

Limit

Western

Blot


Approx 22-23kDa band observed in Mouse Brain and Rat Brain lysates (calculated MW of 22.4kDa according to mouse NP\_032764.1). Recommended concentration: 0.1-0.3µg/ml. An additional very strong band of unknown identity was also consistently observed at 50kDa. This band was successfully blocked by incubation with the immunising peptide. We would appreciate any feedback from people in the field - have any such results been

reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

**Application** 

Pep-ELISA, WB **Type** 







## **GALLERY IMAGES**

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa