

GOAT ANTI-NKD2 ANTIBODY

SKU: EB08176



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	naked cuticle-2 naked cuticle homolog 2 Dvl-binding protein NKD2 naked cuticle homolog 2 (Drosophila) NKD2
Accession ID	NP_149111.1
Blocking Peptide	EBP08176
Immunogen	Peptide with sequence C-HKRYRQKGREGHS, from the internal region of the protein sequence according to NP_149111.1.
Peptide Sequence	C-HKRYRQKGREGHS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Human Gene ID	85409
Mouse Gene ID	72293
Rat Gene ID	308068
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)