

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-NLRP7 ANTIBODY

SKU: EB12896

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

CLR19.4|PYPAF3|NOD12|NALP7|PAN7|HYDM|NLR family, pyrin domain containing 7|NLRP7 Alias

Names

Accession NP_631915.2; NP_996611.2; NP_001120727.1 ID

Blocking

EBP12896 **Peptide**

Peptide with sequence C-QGDIDNFHDDVT, from the internal region of the protein sequence according to Immunogen NP_631915.2; NP_996611.2; NP_001120727.1.

Peptide C-QGDIDNFHDDVT Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human **Gene ID**

199713

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:32000.

Limit

Application Pep-ELISA **Type**