

GOAT ANTI-NLRP7 ANTIBODY

SKU: EB12896



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	CLR19.4 PYPAF3 NOD12 NALP7 PAN7 HYDM NLR family, pyrin domain containing 7 NLRP7
Accession ID	NP_631915.2; NP_996611.2; NP_001120727.1
Blocking Peptide	EBP12896
Immunogen	Peptide with sequence C-QGDIDNFHDDVT, from the internal region of the protein sequence according to NP_631915.2; NP_996611.2; NP_001120727.1.
Peptide Sequence	C-QGDIDNFHDDVT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	199713
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)