

GOAT ANTI-NMDAR2B / GRIN2B ANTIBODY

SKU: EB07294



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	N-methyl-D-aspartate receptor subunit 2B hNR3 NR2B HGNC:4586 glutamate receptor, ionotropic, N-methyl D-aspartate 2B NMDAR2B GRIN2B
Usage Summary	Immunofluorescence: Strong expression of the protein seen in the cytoskeleton/ cytoplasm of HeLa cells. Recommended concentration: 10µg/ml.
Accession ID	NP_000825.1
Blocking Peptide	EBP07294
Immunogen	Peptide with sequence C-KDAHEKDDFHLS, from the internal region of the protein sequence according to NP_000825.1.
Peptide Sequence	C-KDAHEKDDFHLS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human, Rat
Human Gene ID	2904
Mouse Gene ID	14812
Rat Gene ID	24410
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 170kDa band observed in Rat Brain lysates (calculated MW of 166kDa according to NP_000825.1). Recommended concentration: 2-4µg/ml. Primary incubation was 1 hour.
Application Type	Pep-ELISA, WB, IF

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

