

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-NNMT (AA171-182) ANTIBODY

SKU: EB12577



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias NNMT|nicotinamide N-methyltransferase

Names

ID

Accession

NP_006160.1

Blocking

EBP12577 **Peptide**

Peptide with sequence PDLPTYCRALRN, from the internal region of the protein sequence according to **Immunogen**

NP 006160.1.

Peptide Sequence

PDLPTYCRALRN

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human **Species**

Human

4837 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:128000.

Limit

Western Approx 28kDa band observed in lysates of cell line HeLa (calculated MW of 29..6kDa according to

Blot NP_006160.1). Recommended concentration: 0.3-1µg/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES

Telephone: (650) 697-3600



Email: customerservice@vectorlabs.com

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa