

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-NOD1 / CARD4 ANTIBODY

SKU: EB05389



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa 50kDa


37kDa

25kDa

20kDa

15kDa







SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

caspase recruitment domain 4|Nod1|NOD1|caspase recruitment domain family, member 4|CARD4 Names

Accession NP_006083.1 ID

Blocking

EBP05389 **Peptide**

Peptide with sequence EAKVYEDEKRIICF, from the C Terminus of the protein sequence according to **Immunogen**

NP 006083.1.

Peptide Sequence

EAKVYEDEKRIICF

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog, Pig, Cow **Species**

Reactive

Human **Species**

Human

10392 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western Approx 85kDa band observed in Human Duodenum lysates (calculated MW of 107.6kDa according to

Blot NP_006083.1). Recommended concentration: 0.5-1.5µg/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa