


GOAT ANTI-NOD1 / CARD4 ANTIBODY

SKU: EB05389



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	caspase recruitment domain 4 Nod1 NOD1 caspase recruitment domain family, member 4 CARD4
Accession ID	NP_006083.1
Blocking Peptide	EBP05389
Immunogen	Peptide with sequence EAKVYEDEKRIICF, from the C Terminus of the protein sequence according to NP_006083.1.
Peptide Sequence	EAKVYEDEKRIICF
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Reactive Species	Human
Human Gene ID	10392
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 85kDa band observed in Human Duodenum lysates (calculated MW of 107.6kDa according to NP_006083.1). Recommended concentration: 0.5-1.5µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

