

GOAT ANTI-NODAL ANTIBODY

SKU: EB08106



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	nodal, mouse, homolog mouse nodal homolog precursor MGC138230 nodal homolog (mouse) NODAL
Accession ID	NP_060525.2
Blocking Peptide	EBP08106
Immunogen	Peptide with sequence ECPNPVGEEFH, from the internal region of the protein sequence according to NP_060525.2.
Peptide Sequence	ECPNPVGEEFH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Human Gene ID	4838
Mouse Gene ID	18119
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)