

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-NODAL ANTIBODY**

**SKU:** EB08106

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias nodal, mouse, homolog|mouse nodal homolog precursor|MGC138230|nodal homolog (mouse)|NODAL

**Names** 

Accession NP 060525.2 ID

**Blocking** 

EBP08106 **Peptide** 

Immunogen Peptide with sequence ECPNPVGEEFH, from the internal region of the protein sequence according to NP 060525.2.

**Peptide ECPNPVGEEFH** Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** Refrigerated

Instructions

**Predicted** Human, Mouse, Rat **Species** 

Human

4838

Gene ID

Mouse 18119 **Gene ID** 

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:8000.

Limit

**Application** Pep-ELISA

**Type**