

GOAT ANTI-NOS2A ANTIBODY

SKU: EB05258



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	nitric oxide synthase, macrophage NOS, type II HEP-NOS NOS2 INOS NOS nitric oxide synthase 2A (inducible, hepatocytes) NOS2A
Accession ID	NP_000616
Immunogen	Peptide with sequence C-RVAVQPSSLEMSAL, from the C Terminus of the protein sequence according to NP_000616.
Peptide Sequence	C-RVAVQPSSLEMSAL
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	4843
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Jurkat cell lysates at upto 2µg/ml.