




Telephone: (650) 697-3600

GOAT ANTI-NOTCH1 ANTIBODY

SKU: EB07441

e erest



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Alias

Synonym / translocation-associated notch protein TAN-1|notch1|neurogenic locus notch homolog protein 1|Notch (Drosophila) homolog 1 (translocation-associated)|hN1|TAN1|translocation-associated (Drosophila)|Notch

homolog 1|NOTCH1 **Names**

 Immunofluorescence: Strong expression of the protein seen in the nucleus of U2OS cells. Usage

Summary Accession Recommended concentration: 5µg/ml.

ID

NP_060087.2

Blocking Peptide

EBP07441

Immunogen

Peptide with sequence C-DLKARRKKSQDGK, from the internal region of the protein sequence according to NP 060087.2.









Peptide

C-DLKARRKKSQDGK

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method

using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted Species

Human

Reactive

Human

4851

Species Human

Gene ID

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:64000.

Limit

Type

Application Pep-ELISA, IF

GALLERY IMAGES

