Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-NOTCH1A (ZEBRAFISH) ANTIBODY

SKU: EB07958

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias neurogenic locus notch homolog protein|des|deadly seven|n1a|notch homolog 1a|notch1a

Names

Accession NP 571516.1

Blocking

EBP07958 **Peptide**

Peptide with sequence C-QNMAPMQHQNISHH, from the internal region of the protein sequence according to Immunogen

NP 571516.1.

Peptide C-QNMAPMQHQNISHH Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Zebrafish

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection

Antibody detection limit dilution 1:64000.

Limit

Western Not yet tested - our routine western blotting protocol does not allow for the detection of proteins >250kDa **Blot** (calculated size of approx. 262kDa according to NP_571516.1) and the antibody is an aspiring product.

Application

Pep-ELISA **Type**