

GOAT ANTI-NOTCH1A (ZEBRAFISH) ANTIBODY

SKU: EB07958



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	neurogenic locus notch homolog protein des deadly seven n1a notch homolog 1a notch1a
Accession ID	NP_571516.1
Blocking Peptide	EBP07958
Immunogen	Peptide with sequence C-QNMAPMQHQNISHH, from the internal region of the protein sequence according to NP_571516.1 .
Peptide Sequence	C-QNMAPMQHQNISHH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Zebrafish
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Not yet tested - our routine western blotting protocol does not allow for the detection of proteins >250kDa (calculated size of approx. 262kDa according to NP_571516.1) and the antibody is an aspiring product.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)