

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-NOTCH1B (ZEBRAFISH) ANTIBODY

SKU: EB07959

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias notch receptor protein 1b|etlD309871.5|wu:fc55e11|fc55e11|notch homolog 1b|notch1b

Names

Accession NP 571377.2 ID

Blocking

EBP07959 **Peptide**

Peptide with sequence C-HLGIGNNQDLSHIQ, from the internal region of the protein sequence according to Immunogen

NP 571377.2.

Peptide

C-HLGIGNNQDLSHIQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Zebrafish

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection

Western

Antibody detection limit dilution 1:64000.

Limit

Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We

Blot would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Application

Pep-ELISA

Type