

GOAT ANTI-NPC1L1 ANTIBODY

SKU: EB09948



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	NPC1L1 NPC1L1 NPC1 (Niemann-Pick disease, type C1, gene)-like 1 Niemann-Pick C1-like protein 1
Accession ID	NP_037521.2; NP_001095118.1
Blocking Peptide	EBP09948
Immunogen	Peptide with sequence CNNFSFTQKIYAT, from the internal region of the protein sequence according to NP_037521.2; NP_001095118.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_037521.2; NP_001095118.1).
Peptide Sequence	CNNFSFTQKIYAT
Shipping Instructions	Refrigerated
Predicted Species	Human, Cow
Human Gene ID	29881
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Preliminary experiments in lysates of cell lines HeLa, Jurkat and K562 gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA