

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-NPC1L1 (MOUSE) (AA56-69) ANTIBODY

SKU: EB10837

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

OTTMUSP0000000725| Niemann-Pick disease, type C1-like 1| Niemann-Pick disease, type C1| Niemann-Pick C1-

like protein 1| Gm243| 9130221N23Rik| RP23-198N14.2| NPC1-like 1|Npc1|1

Names Accession

Alias

NP 997125.2

Blocking

EBP10837 **Peptide**

Peptide with sequence CLSNTPARHVTGDH, from the N Terminus of the protein sequence according to Immunogen

NP 997125.2.

Peptide

CLSNTPARHVTGDH Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Mouse, Rat

Mouse Gene ID

237636

Rat Gene ID 432367

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection

Antibody detection limit dilution 1:128000.

Limit

Western Blot

Preliminary experiments in Mouse and Rat Liver lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been

reported with other antibodies/lysates?

Application

Type

Pep-ELISA





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600