

GOAT ANTI-NPHP1 ANTIBODY

SKU: EB05453

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	SLSN1 NPHP1 nephronophthisis 1 (juvenile) NPHP1
Names	
Accession ID	NP_000263
Immunogen	Peptide with sequence C-TYDFLGEMRKNV, from the C Terminus of the protein sequence according to NP_000263.
Peptide Sequence	C-TYDFLGEMRKNV
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Mouse
Human Gene ID	4867
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 75-80kDa band seen in Mouse Kidney lysate. Recommended for use at 1-3µg/ml.