

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-NPHP1 ANTIBODY**

**SKU:** EB05453

×

## **SPECIFICATIONS**

 $\textbf{Formulation} \ \text{Supplied at 0.5 mg/ml in Tris saline, 0.02\% sodium azide, pH7.3 with 0.5\% bovine serum albumin.}$ 

 $\textbf{Unit Size} \qquad 100 \; \mu \text{g}$ 

Storage

Instructions

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias SLSN1|NPH1|nephronophthisis 1 (juvenile)|NPHP1

Names

Accession NP\_000263

Peptide with sequence C-TYDFLGEMRKNAV, from the C Terminus of the protein sequence according to

Immunogen NP 000263.

**Peptide** 

Sequence C-TYDFLGEMRKNAV

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

**Method** using the immunizing peptide.

Shipping

**Instructions** Refrigerated

Predicted

Species Human

Reactive

Species

Mouse

Human

Gene ID 4867

Product

Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

Western

Western blot: Approx 75-80kDa band seen in Mouse Kidney lysate. Recommended for use at 1-3μg/ml.