Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-NRAP ANTIBODY

SKU: EB10399

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias NRAP|nebulin-related anchoring protein

Names

Accession NP 006166.3; NP 932326.2 ID

Blocking

EBP10399 **Peptide**

Peptide with sequence C-DFAKSRSQFH, from the internal region of the protein sequence according to Immunogen

NP 006166.3; NP 932326.2.

Product This antibody is expected to recognize both reported isoforms (NP_006166.3; NP_932326.2). Comments

Peptide

C-DFAKSRSQFH Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Refrigerated

Instructions

Predicted Human, Mouse, Rat, Dog, Cow, Pig Species

Human

4892

Gene ID Mouse 18175 **Gene ID**

Rat Gene ID 307982

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Application Pep-ELISA Type





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600