

GOAT ANTI-NRAP ANTIBODY

SKU: EB10399



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	NRAP nebulin-related anchoring protein
Accession ID	NP_006166.3; NP_932326.2
Blocking Peptide	EBP10399
Immunogen	Peptide with sequence C-DFAKSRSQFH, from the internal region of the protein sequence according to NP_006166.3; NP_932326.2.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_006166.3; NP_932326.2).
Peptide Sequence	C-DFAKSRSQFH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow, Pig
Human Gene ID	4892
Mouse Gene ID	18175
Rat Gene ID	307982
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)