

GOAT ANTI-NRARP ANTIBODY

SKU: EB12330



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. |
| Unit Size | 100 µg |
| Storage Instructions | Aliquot and store at -20°C. Minimize freezing and thawing. |
| Synonym / Alias Names | notch-regulated ankyrin repeat-containing protein NOTCH-regulated ankyrin repeat protein NRARP |
| Accession ID | NP_001004354.1 |
| Blocking Peptide | EBP12330 |
| Immunogen | Peptide with sequence SQAELSTCSAPQTQR, from the N Terminus of the protein sequence according to NP_001004354.1. |
| Peptide Sequence | SQAELSTCSAPQTQR |
| Purification Method | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| Shipping Instructions | Refrigerated |
| Predicted Species | Human, Mouse, Rat, Dog, Pig, Cow |
| Human Gene ID | 441478 |
| Mouse Gene ID | 67122 |
| Rat Gene ID | 499745 |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png |
| ELISA Detection Limit | Antibody detection limit dilution 1:1000. |
| Application Type | Pep-ELISA |

DOCUMENTS

- [Data Sheet](#)