Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-NRARP ANTIBODY**

**SKU:** EB12330

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias notch-regulated ankyrin repeat-containing protein|NOTCH-regulated ankyrin repeat protein|NRARP

**Names** 

Accession NP 001004354.1

ID

**Blocking** 

EBP12330 **Peptide** 

Peptide with sequence SQAELSTCSAPQTQR, from the N Terminus of the protein sequence according to Immunogen

NP 001004354.1.

**Peptide** 

**SQAELSTCSAPQTQR** Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

**Species** 

Human, Mouse, Rat, Dog, Pig, Cow

Human

441478 Gene ID

Mouse **Gene ID** 

67122

**Rat Gene ID** 499745

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:1000.

Limit

**Application** Pep-ELISA

**Type**