



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-NSP3 ANTIBODY

SKU: EB05105

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

likely ortholog of mouse Cas/HEF1-associated signal transducer|novel SH2-containing protein 3|SH2 domain-Alias

containing 3C|SH2 domain containing 3C|NSP3|CHAT|SH2D3C

Names

Accession NP 005480

Blocking

EBP05105 **Peptide**

Immunogen Peptide with sequence C-SHKLEPAVRSSEL, from the C Terminus of the protein sequence according to NP 005480.

Peptide

C-SHKLEPAVRSSEL Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

10044

Gene ID

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

Western

Western blot: Band observed in Jurkat at approx 50kDa with 1 hour incubation of primary. Blot