

GOAT ANTI-NSP3 ANTIBODY

SKU: EB05105



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	likely ortholog of mouse Cas/HEF1-associated signal transducer[novel SH2-containing protein 3 SH2 domain-containing 3C SH2 domain containing 3C NSP3 CHAT SH2D3C
Accession ID	NP_005480
Blocking Peptide	EBP05105
Immunogen	Peptide with sequence C-SHKLEPAVRSSSEL, from the C Terminus of the protein sequence according to NP_005480.
Peptide Sequence	C-SHKLEPAVRSSSEL
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	10044
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: Band observed in Jurkat at approx 50kDa with 1 hour incubation of primary.