



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-NTE ANTIBODY

SKU: EB09135

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias neuropathy target esterase|sws|SPG39|NTEMND|NTE|patatin-like phospholipase domain containing 6|PNPLA6

Names

Accession NP 006693.3 ID

Blocking

EBP09135 **Peptide**

Peptide with sequence C-YQGRQTSSIFEAAKQ, from the internal region of the protein sequence according to Immunogen

NP 006693.3.

Peptide

C-YQGRQTSSIFEAAKQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse

Human

10908 Gene ID

Mouse

50767

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:64000.

Limit

Application Pep-ELISA **Type**