

GOAT ANTI-NTE ANTIBODY

SKU: EB09135



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	neuropathy target esterase sws SPG39 NTEMND NTE patatin-like phospholipase domain containing 6 PNPLA6
Accession ID	NP_006693.3
Blocking Peptide	EBP09135
Immunogen	Peptide with sequence C-YQGRQTSSIFEAAKQ, from the internal region of the protein sequence according to NP_006693.3.
Peptide Sequence	C-YQGRQTSSIFEAAKQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	10908
Mouse Gene ID	50767
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)