



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-OAS2 ANTIBODY

SKU: EB12372



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

OAS2|2'-5'-oligoadenylate synthetase 2, 69/71kDa|(2'-5')oligo(A) synthetase 2|2'-5'-oligoadenylate synthase

Alias 2|2-5A synthase 2|p69 OAS / p71 OAS Names

Accession

NP_058197.2; NP_002526.2; NP_001027903.1 ID

Blocking

EBP12372 **Peptide**

Peptide with sequence C-RKTVLRGNSDGT, from the internal region (near N terminus) of the protein sequence **Immunogen**

according to NP_058197.2; NP_002526.2; NP_001027903.1.

Product Comments

This antibody is expected to recognize all reported isoforms (NP 058197.2; NP 002526.2; NP 001027903.1).

Peptide

Sequence

C-RKTVLRGNSDGT

4939

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human **Species**

Human

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:128000.

Limit

Western Approx 85 +19kDa bands observed in lysates of cell line Jurkat (calculated MW of 82.4kDa according to NP_058197.2 and of 19.5kDa according to NP_001027903.1). Recommended concentration: 1-3μg/ml. Blot

Application

Pep-ELISA, WB Type

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa